

CASTELLINI

EN

TREATMENT
C7 C8 C9
Stools



MAKING
SPACE
FOR
EXCELLENCE

C7 C8 C9 Stools

COMFORT AND WELLBEING

Designed to ensure correct posture and wellbeing in every working position, the latest dentist and assistant stools are designed to fit effortlessly into any dental practice.

Three versions of comfort.




C7 STOOL

The C7 is an optimised stool that lets dentists perform treatment perfectly while enjoying correct body support. Adjustment mechanisms ensure both seat and back are always regulated to suit the operator's build, while the compact design allows full access to all parts of the dental unit. The lining material and the manufacturing technique make it easy to sanitise the surfaces.



SEAT

The seat is specially shaped to aid natural positioning of the legs, ensuring a stable, comfortable posture at all times. Soft, comfortable foam padding contributes to the dentist's wellbeing, even during prolonged treatment sessions.



BACKREST

The small backrest reduces the contact surface area and, therefore, the sensation of heat. Simultaneously, its shape provides perfect lumbar support during treatment or when talking with the patient.



SEAT ADJUSTMENT

Seat height can be adjusted as desired to ensure proper leg support and spinal posture at all times. An optional accessory can add a further 10 cm of height.



LUMBAR SUPPORT

The backrest has 9 different height settings, providing proper support for all dentists whatever their build.



BACKREST ADJUSTMENT

The backrest can be tilted as preferred by the user, who can always 'model' the stool to suit their build.

C9 STOOL

Innovative saddle-like seat design lets dentists approach the patient without having to flex their back, reducing the risk of fatigue. Maximum comfort with the smallest possible footprint.



IDEAL POSTURE

Saddle-like shaping angles the legs correctly, consequently improving back posture and making it easier to get close to the treatment zone. Excellent ergonomics ensures a feeling of wellbeing and energy, minimising any risk of musculoskeletal pain.



DYNAMIC SEATING

A tilt mechanism ensures the 'saddle' supports the dentist's body movements, providing dynamic seating to ensure a proper, active posture in every treatment scenario. The angle of movement is easily adjusted and can be locked in place.



UPHOLSTERY

Upholstery is available in two versions: seamless upholstery to make sanitising tasks easier and Memory Foam padding (optional) to provide very high comfort levels.

C8 STOOL

The C8 stool is designed to provide maximum mobility, quick access to all instruments and optimal arm support.



SYMMETRICAL SEAT

The seat is circular, making it comfortable even when operators rotate their body. Seat and backrest are actually free to move so the operator can move through 360°. The armrest framework is designed to keep the forearm stable in any position.



COMFORTABLE, PRACTICAL ARMREST

The armrest is horizontal, without edges and lined with non-slip material: assistants can therefore keep their arm in a comfortable, stable position, even during long treatment sessions. During the assembly

stage, users can opt for a centred armrest position, or off-centred to the left or right according to the needs of the practice.

TECHNICAL SPECIFICATIONS


C7 Stool



C8 Stool



C9 Stool



Measurements in drawings are in millimetres

	Footrest ring	10 cm extension	Balance system	Standard upholstery	Memory Foam padding
C7	○	○	-	●	-
C8	○	○	-	●	-
C9	○	-	●	●	○

● = supplied as standard ○ = optional - = not available



EXTENSION

An optional 10 cm extension is available for all models, ensuring the seat can be adapted to all user builds.



FOOTREST RING

An optional footrest ring providing comfortable, stable support is available on every model. A metallic finish ensures easy cleaning.



STANDARD

- 102 Atlantic Blue
- 106 Mediterranean Blue
- 113 Pacific Blue
- 136 Indian Blue
- 132 Blueberry Violet
- 134 Japanese Wisteria
- 135 Venetian Red
- 115 Scottish Salmon
- 103 Nevada Yellow
- 123 Polynesian Green

- 101 Caribbean Green
- 137 Bright Silver
- 121 Anthracite Grey
- 130 Graphite Black
- 142 Anatolian Hazelnut
- 141 Papyrus Beige
- 143 Arabian Gold
- 140 Brazilian Brown
- 144 Ruby Red

MEMORY FOAM

- 198 Atlantic Blue
- 196 Mediterranean Blue
- 183 Pacific Blue
- 186 Indian Blue
- 192 Blueberry Violet
- 184 Japanese Wisteria
- 194 Venetian Red
- 195 Scottish Salmon
- 182 Nevada Yellow
- 193 Polynesian Green

- 197 Caribbean Green
- 187 Bright Silver
- 199 Anthracite Grey
- 180 Graphite Black
- 172 Anatolian Hazelnut
- 171 Papyrus Beige
- 173 Arabian Gold
- 170 Brazilian Brown
- 174 Ruby Red

CASTELLINI

PASSION FOR DENTISTRY

Since 1935

CSTOGB261500 02/2026

Due to our policy of constant technological upgrading, technical specifications may be subject to change without prior notice. According to the standards in force, in extra-EU areas the availability and specifications of some products and/or characteristics may vary. Please contact your local distributor for further information.



**Bu Medical Equipment
Sede Legale Ed Amministrativa
Headquarters**

Cefla s.c.
Via Selice Provinciale, 23/a - 40026
Imola - Bo (Italy)
tel. +39 0542 653111
fax +39 0542 653344

castellini.com

**Stabilimento
Plant**

Via Bicocca, 14/c - 40026
Imola - Bo (Italy)
tel. +39 0542 653441
fax +39 0542 653601